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# (EN) Table of contents

<b>1. Foreword</b> .....	<b>17</b>
1.1 Information about the manual.....	17
1.2 Notes on the operating instructions .....	17
1.3 Representation of safety instructions.....	17
1.3.1 Safety instructions in the operating instructions .....	17
<b>2. Safety</b> .....	<b>17</b>
2.1 General .....	17
2.2 General information regarding safety and warning instructions .....	17
<b>3. Residual risks and foreseeable misuse</b> .....	<b>18</b>
<b>4. Description</b> .....	<b>18</b>
4.1 Description of the unit .....	18
4.2 Scope of delivery.....	18
4.3 Correct use.....	19
4.4 Technical data.....	19
<b>5. Preparation and commissioning</b> .....	<b>19</b>
5.1 Preparation .....	19
5.1.1 Positioning the AirOne .....	19
5.1.2 Nebulising chamber and nebuliser cups .....	19
5.1.3 Preparing the horse for inhalation .....	19
5.2 Commissioning .....	19
5.2.1 Adjusting the fresh air supply on the breathing mask .....	19
5.2.2 Meaning of the LED displays .....	20
5.2.2.1 Temperature LED lights up .....	20
5.2.2.2 Water level LED lights up.....	20
<b>6. Maintenance and care</b> .....	<b>20</b>
6.1 Careful cleaning and disinfection.....	20
6.2 Important information about the contact liquid .....	20
6.3 Further care instructions .....	20
<b>7. Troubleshooting</b> .....	<b>21</b>
<b>8. Warranty conditions</b> .....	<b>21</b>
<b>9. CE-/UKCA-conformity declaration</b> .....	<b>21</b>
<b>10. Electrical waste</b> .....	<b>21</b>
<b>11. Spare parts and accessories</b> .....	<b>22</b>

# (EN) Read and pay attention to these operating instructions before assembly and initial commissioning!

## Keep these operating instructions for future reference!

### 1. Foreword

#### 1.1 Information about the manual

We reserve the right to make changes due to further technical development compared to the data and illustrations mentioned in this manual. Reprints, translations and reproductions in any form, including excerpts, require the written consent of the manufacturer.

This manual is part of the scope of delivery.

- It must be kept close at hand and remains with the unit even if it is sold.
- This manual is not subject to a revision service. The respective current status can be obtained from the specialist trader or directly from the manufacturer.
- It has a modular structure and refers exclusively to the named product.
- These operating instructions describe the operation, use and maintenance of the unit.

The safety instructions are intended to prevent personal injury, damage to animals and the unit. All operators are required to read and observe these safety instructions at all times. The manufacturer accepts no liability or warranty for damage and malfunctions resulting from non-observance of the operating instructions!

#### 1.2 Notes on the operating instructions

Version/revision: 1.0

Creation date: 2021-11-25

Last change: 2022-04-07

Language of the original operating instructions: German

#### 1.3 Representation of safety instructions

##### 1.3.1 Safety instructions in the operating instructions

#### **Danger!**

Imminent danger which leads to serious physical injury or death.

#### **Caution!**

Risk of physical injury to people and animals and material damage if used incorrectly!

#### **Important!**

Potentially harmful situation that could damage the product or an object nearby.

#### **Note**

This pictogram identifies information that contributes to a better understanding of the work processes.

### 2. Safety

#### 2.1 General

- Read the following instructions carefully before commissioning and using the unit.
- Anyone tasked with commissioning, operation or maintenance must have read and understood the operating instructions and in particular the chapter on safety.
- The unit must only be operated, maintained and repaired by trained and authorised personnel.
- Failure to comply with the individually listed points in the operating instructions can result in physical injury or material damage for which the manufacturer assumes no liability.

#### 2.2 General information regarding safety and warning instructions

#### **Danger!**

**Imminent danger which leads to serious physical injury or death.**

- Only open the unit when unplugged!
- Work on the electrical cables may only be carried out by a qualified electrician.
- In the event of a fault, the mains plug must be disconnected first.
- The unit must be operated from a 230 V socket outlet.
- Only extension cables that have a suitable safety class and splash protection may be used to connect the unit

#### **Caution!**

**Risk of physical injury to people and animals and material damage if used incorrectly!**

- Children and infirm individuals must only use the device under supervision.
- Note the safety and maintenance instructions.

- The unit must be used exclusively in accordance with intended use and in a technically safe condition! The operational safety of the unit is only guaranteed if it is used for its intended purpose. Failure to comply with the individually listed points in the operating instructions can result in physical injury or material damage for which the manufacturer assumes no liability.
- Keep all safety and danger notices on/at the unit complete and always in a legible condition.
- The maintenance intervals specified in these operating instructions must be observed.
- Only use genuine spare parts.
- If the **mains cable** for this device is damaged, then to avoid risk the cable must be replaced by the manufacturer, its customer service team, or a similarly qualified person.

### 3. Residual risks and foreseeable misuse

#### Risks for persons:

under the following circumstances, the safety of persons and material objects can be endangered:

- Troubleshooting, adjustment, cleaning, repair and maintenance work contrary to the instructions in the operating instructions.
- If untrained personnel are working on electronic parts.
- Failure to comply with the safety instructions in the operating instructions.
- Unauthorised modifications to the unit.
- Use of non-original manufacturer spare parts.

#### The following points summarise further measures aimed at protecting against personal and material risks:

- Mount the unit in a place that is well ventilated and free from electromagnetic interference.

#### Risk for animals:

- The animals are not yet accustomed to using the inhaler (risk of panic reaction).
- The animals are unsupervised during use.
- The animals are not tethered correctly to the site.

## 4. Description

### 4.1 Description of the unit

#### Ultrasound Inhaler for Horses

The devices are designed to split (nebulise) inhalants by ultrasound into components small enough for them to be optimally absorbed by the horse through inhalation. The choice of inhalant is the operator's responsibility.

#### Caution!

The unit must be used exclusively in accordance with intended use and in a technically safe condition! The operational safety of the unit is only guaranteed if it is used for its intended purpose. Failure to comply with the individually listed points in the operating instructions can result in physical injury or material damage for which the manufacturer assumes no liability.

### 4.2 Scope of delivery

#### The AirOne ultrasound inhaler is made up of the following components:

1. Ultrasound unit incl. mains connection cable
2. Nebulising chamber (fixed component of the nebuliser)
3. Removable nebulising chamber cover with connecting pieces (3.1)
4. Breathing mask with connecting pieces (4.1), adjustable fresh air supply (4.2), rubber seal (4.3), adjustable head strap (4.4) and one inhalation (4.5) and two removable exhalation valves (4.6)
5. Breathing tube
6. Wall bracket with 4 screws, 4 dowels and installation instructions
7. Transport box with cover (not shown)

#### The accessories also include:

a measuring cup, 30 nebuliser cups, 1 litre of contact liquid.

### 4.3 Correct use

This unit must only be used to generate ultrasonic mist with the original contact liquid. The unit must be used for its intended purpose with an original Hippomed mask of the correct size exclusively on horses. The use of nebuliser add-ons is the sole responsibility of the operator. The regulations set out in these operating instructions must be complied with. Any other use of the unit, e.g. in any other area than those described, shall be deemed incorrect use. The manufacturer accepts no liability for any resulting physical and/or material damage.

### 4.4 Technical data

Item number / type designation: AirOne ultrasound nebuliser  
Nebuliser output: Average: 6.7 ml/min, Maximum: 15 ml/min  
Particle size: 0.47 - 6 µm  
Mains supply: AC 230 V 50 Hz  
Current consumption: 300 mA  
Operating mode: Short-term operation KB, max. 30 minutes, then manually set aside to cool down  
Protection class: II  
Protection type: IP 54  
Weight of base unit with cable: 1560 g

#### Ambient conditions

during operation: 0°C to 40°C and 700 to 1060 hPa

## 5. Preparation and commissioning

Congratulations for choosing a quality product from Hippomed! The AirOne features animal-friendly nebulisation and technical sophistication. The AirOne is a dual-chamber, ultrasonic cold nebulisation system. The AirOne was developed especially for horses and has been immensely popular for over 20 years with veterinary surgeons, veterinary practices, breeders and stable owners, as well as private horse lovers.

### 5.1 Preparation

#### 5.1.1 Positioning the AirOne

Please attach the AirOne's wall bracket as shown in the installation instructions supplied. Should you not yet have decided on a permanent location for the unit, place the AirOne on a table, trolley or similar, but never on the floor. The raised position of around 120 to 150 cm guarantees a high degree of effectiveness for the nebulisation process.

#### 5.1.2 Nebulising chamber and nebuliser cups

Grip the nebulising chamber cover at the tube connection piece and pull it upwards with a slight tilting movement. The locking mechanism releases and the nebulising chamber cover can be removed and placed to one side. Pour in **60 ml of contact fluid** into the nebulising chamber using the measuring cup. Please shake the contact fluid well before doing so. Then place the nebuliser cup in the nebulising chamber and add the medium to be nebulised to the cup up to a **maximum fill level of 50 ml**. Replace the nebulising chamber cover and push it down firmly until it closes flush with the upper edge of the housing.

**Important:** the AirOne delivers its best performance with a **MAXIMUM** fill level of **50 ml** in the nebuliser cup. **You should therefore avoid exceeding this amount.**

The duration of inhalation is determined by the properties of the preparation to be nebulised. Please do not nebulise any **oil-containing** substances. The generally recommended duration for inhalation is around 15 - 20 minutes. For technical reasons, the nebulisation output diminishes towards the end of a treatment session. A residual volume of around 5 ml may remain in the nebuliser cup. A comparatively long nebulising session is needed to nebulise this remaining amount. We therefore recommend ending the treatment session.

**Also, always make sure** that the white clamping ring in the attachment piece of the nebulising chamber cover is located as far as possible at the front of the opening towards the tube and that the slot in the clamping ring is facing upwards. This prevents too much of the bubbling medium from travelling into the tube and accumulating there. Condensation may occur in the breathing tube during treatment. Occasional raising of the tube can cause a backflow of the condensation water into the nebuliser cup.

After you have connected the supplied tube to the connection piece of the breathing mask and the nebulising chamber cover, the AirOne is ready for operation once you have plugged the mains plug into the earthed socket.

#### 5.1.3 Preparing the horse for inhalation

For safety reasons, the horse being treated must wear a halter and be tethered on both sides with rubber ropes (www.kerbl.de, Art No. 3210370 - 3210372). Check the correct fit of the rubber seal on the breathing mask. The arrow on the rubber seal must point towards the front of the mask and down towards the connection piece. Push the mask with the rubber ring over the horse's nostrils and fasten it securely with the head strap.

### Caution!

**The fresh air supply must be fully open to avoid possible panic reactions by the horse.**

### 5.2 Commissioning

Switch the AirOne on by pressing the ON/OFF button on the control panel. The function is signalled by a green operation indicator lighting up. Inhalation of the mist is best checked on the top of the nebulising chamber cover.

#### 5.2.1 Adjusting the fresh air supply on the breathing mask

Stepless closure of the fresh air supply on the breathing mask means that the nebulisation quantity and thus the efficiency can be varied. However, this should only be done very carefully and with simultaneous observation of the horse's behaviour. The animal should show calm, even breaths. If the fresh air supply on the breathing mask is closed completely too soon, the horse may panic. The horse must first learn inhalation (breathing) and therefore the inhalation times vary depending on the horse and its state of health. The inhalation time is 10 to max. 20 minutes.

## 5.2.2 Meaning of the LED displays

### 5.2.2.1 Temperature LED lights up

When the "Temperature" indicator on the AirOne's control panel (membrane keypad) lights up, the unit has switched off due to overheating of the contact liquid. Disconnect the mains plug and allow the contact liquid to cool down for 10 minutes.

**i** The warmer the contact medium, the sooner the unit switches off. The same also applies to extremely high outside temperatures, in which case you should use a cooled contact medium.

#### Caution!

To avoid a hazard due to an unintentional reset of the protective temperature limiter, the unit must not be supplied by an external switching device such as a timer or be connected to a circuit that is regularly switched on and off by a device.

### 5.2.2.2 Water level LED lights up

When the "Water level indicator" lights up on the AirOne's control panel (membrane keypad), there is normally too little contact liquid in the nebulising chamber. Switch off the unit, remove the nebulising chamber cover and nebuliser cup and replace the contact liquid. After switching on, treatment with the AirOne can be continued. If the control indicator continues to light up, this could also be due to deposits on the sensors in the nebulising chamber. Carefully remove the deposits with a soft brush and a little vinegar-based cleaner. If the unit still does not work after this, there is probably a technical defect and you must send it to the service centre for inspection.

#### Caution!

**Unplug the unit from the mains before any maintenance/cleaning work.**

## 6. Maintenance and care

### 6.1 Careful cleaning and disinfection

To ensure a sustained high level of treatment success, it is essential that the AirOne is carefully cleaned and disinfected after each use. Commercially available washing-up liquid is suitable for external cleaning. For disinfection, please use the disinfectant included with the AirOne. Follow the instructions for use on the disinfectant label. Remove the exhalation valves on the breathing mask by turning them a quarter turn to the left (bayonet catch) and remove the nebulising chamber cover on the nebuliser. Carefully spray all components, namely the mask, valves and nebulising chamber cover. Hang up the breathing tube and spray into the tube. Please observe the information on exposure time on the disinfectant's instructions for use. The nebuliser cup can be used again after cleaning. However, it is recommended that the cup is replaced after 3 treatments.

After the exposure time, rinse all parts with fresh water and leave to dry in a suitable place. Then reassemble all the parts and store them in the transport box provided. The AirOne is thus perfectly prepared for the next use.

### 6.2 Important information about the contact liquid

#### Caution!

**Long-term trouble-free operation is only guaranteed if the contact fluid supplied is used.** If other agents such as distilled or tap water are used as a contact medium, the oscillator unit, the heart of the unit, will be **permanently damaged** and the unit will then have to be repaired and all warranty claims will be invalidated. The contact agent is produced in the laboratory with a conductivity specially adapted to the AirOne. If this is not or no longer correct, the water level indicator lights up.

**i** The contact liquid is only renewed when the water level indicator lights up or there is visual contamination. If the AirOne is not used for more than a week, please remove the contact fluid and only refill the AirOne when it is used again.

### 6.3 Further care instructions

#### Caution!

**Do not carry out any descaling activities. These permanently damage the surface at the bottom of the nebulising chamber (ultrasonic oscillator / transducer) and it will need to be replaced. Never touch or clean the oscillator with a hard object. Do not use acidic cleaning agents.** For optimal care of the nebulising chamber, take a soft brush, pour water without additives into the nebulising chamber and use the brush to swirl up any possible residue and then rinse it out. You can wipe the housing with a damp cloth and then rub it dry.

## Any questions?

We'll be happy to advise you - simply contact us on:

T: +49 (0)8086 933 327

T: +49 (0)5204 92970

E: inhalatoren@kerbl.com

For chargeable or warranty repairs, please download the Hippomed repair application and fill it out in full:

<https://www.kerbl.com/de/service/downloads>

## 7. Troubleshooting

Error	Remedy
Temperature LED lights up	Chapter 5.2.2.1
Water level LED lights up	Chapter 5.2.2.2
Horse is not inhaling properly	Check the inhalation valves for contamination Check the tube connection between the unit and mask Check the air supply setting on the mask
No inhalation mist visible	Check the content of the nebuliser cup  Unit not switched on (switch LED must light up green)  Oscillator contaminated - see Chapter 6.3  Oscillator faulty - contact the service centre
Unit does not switch on	No power - check the power supply  Switch faulty - contact the service centre

## 8. Warranty conditions

Albert Kerbl GmbH grants you a 2-year warranty on the ultrasound base unit in accordance with our general warranty conditions. These general warranty conditions can be found at [www.kerbl.com/guarantee](http://www.kerbl.com/guarantee). We will also be happy to send you the warranty conditions by post. Your statutory warranty rights remain unaffected by our guarantee.

This warranty begins on the date of purchase. The warranty does not cover:

1. Damage caused by natural wear and tear (wear parts such as the oscillator / transducer), excessive strain and incorrect use (non-compliance with operating instructions).
2. Damage caused by force majeure or sustained during transport.
3. Damage to rubber parts, e.g. on access lines.
4. Damage to the mask including its attachments (valves) and the tube.

**Warranty repairs are only carried out by the service centre (DE/AT) or an authorised contracted workshop! The manufacturer reserves the right to reject any warranty claims in the event of non-compliance.**

The unit must for this purpose be delivered to the specified service address, at the risk of the unit's owner, complete, packaged and with sufficient postage, including the original invoice.

Please provide as much information as possible on the nature of the problem (see form).

**Costs incurred as a result of unjustified defect claims will be invoiced to the unit's owner.**

### 9. CE-/UKCA-conformity declaration

Albert KERBL GmbH hereby declares that the product / device described in these instructions complies with the fundamental requirements and other relevant stipulations and regulations. The CE-/UKCA-mark confirms compliance with the Directives of the European Union or the relevant UK legislation.

### 10. Electrical waste

It is up to the operator to dispose of the device responsibly after its useful life in accordance with local regulations. Consult the relevant regulations for your country. The device must not be disposed of in household waste. Under the terms of the EU Directive on the disposal of old electrical and electronic equipment, the unit can be taken to local authority-run refuse sites or recycling plants free of charge, or returned to dealers who offer a returns service. Correct disposal protects the environment and prevents any possible harmful impacts on humans and the environment.

## 11. Spare parts and accessories

Item number	Article designation	Illustration
3211802	Complete inhaler with warmblood mask and accessories in transport box	—————
3211800	Inhaler with accessories in transport box, without mask	• • • • • • • • • •
3211800-0	Inhaler without accessories	1,2,3,3.1
3211821	Wall bracket incl. 4 screws and 4 dowels with installation instructions	6
3211826	Transport box with cover	7
3211807	Complete warmblood mask with valves, seal and head straps	4-4.6
3211809	Complete coldblood mask with valves, seal and head straps	4-4.6
3211810	Complete Shetty mask with valves, seal and head straps	4-4.6
3211827	Coldblood mask without accessories	4
3211829	Shetty mask without accessories	4
3211808	Warmblood mask without accessories	4
3211832	Seal for warmblood mask	4.3
3211833	Seal for Shetty mask	4.3
3211859	Seal for coldblood mask	4.3
3211830	Inhalation valve for breathing mask	4.5
3211831	Exhalation valve for breathing mask	4.6
3211825	Head strap for breathing mask	4.4
3211823	Breathing tube 170 length	5
3211839	Nebulising chamber cover	3
3211841	Snap ring for nebulising chamber cover	8
3211818	Nebuliser cup, pack of 30	10
3211822	Measuring cup 250 ml	9
3211819	Contact fluid 1 litre	11
299310	Disinfectant, 1 litre(Please ask for availability outside DE/AT!)	12

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### **Kerbl Austria Handels GmbH**
### **Kerbl UK Ltd**