

## GLASS MEASURING CYLINDER



### DESCRIPTION

Stable measuring cylinder made of high-quality borosilicate glass; available in various sizes. The scale on the side is easy to recognise with its blue colour and eases liquid measurement

### PRODUCT DETAILS

- Glass measuring cylinder
- Material: 3.3 borosilikate glass
- Easy-to-read blue scale
- With spout and hexagonal glass base
- Good stability
- Available in the following sizes: 100 ml, 250 ml and 500 ml
- Delivered quantity: 1 piece




### Overview

Article Number:  
**T135854-T135856**

Category:  
**Lab Equipment**

Teqler is a trademark of NetMed S.à.r.l.

 NetMed S.à.r.l.  
8, Am Scheerleck, 6868 Wecker, Luxembourg

 +352 267149 09  +352 267149 19  [info@netmed.lu](mailto:info@netmed.lu)

Social Media

