

FIELD BED



DESCRIPTION

Collapsible camp bed in the colour dark green with robust cover and 6 stable feet. Ideal for outdoor use – e.g. medical service or disaster control.

PRODUCT DETAILS

- Camp bed, dark green
- Only 8,5kg heavy
- Foldable to 90 x 12 x 20cm
- Loadable up to 159kg
- Robust cover
- 6 feet
- CE mark
- FDA certified
- Size in assembled condition: 194 x 66 x 41cm
- Transport bag included in scope of delivery
- 1 piece

Overview

Article Number:
T138244



Category:
Emergency Equipment

EAN:
4260306772032



Teqler is a trademark of NetMed S.à.r.l.

 NetMed S.à.r.l.
8, Am Scheerleck, 6868 Wecker, Luxembourg

 +352 267149 09  +352 267149 19  info@netmed.lu

Social Media

