

## TEQLER PRECISION DOSING SYRINGE



### DESCRIPTION

This fine dosage syringe is specially designed for the administration of particularly small amounts of medication. The fine scale ensures precise winding, the economy spike prevents the loss of medication.

### PRODUCT DETAILS

- 3-part disposable syringe with 1ml volume
- with economy spike
- Scaling: 0.01ml
- Luer connector
- Smooth running
- Sealed individually sterile
- 1 box of 100 pieces

### Overview

Article Number:  
**T138322**




Category:  
**Infusion | Injection | Syringes**

EAN:  
**4260626263746**



Teqler is a trademark of NetMed S.à.r.l.

 NetMed S.à.r.l.  
8, Am Scheerleck, 6868 Wecker, Luxembourg

 +352 267149 09  +352 267149 19  [info@netmed.lu](mailto:info@netmed.lu)

Social Media

