

enbio 

# Magic Filter

Quick Start

Ver.3 04.12.2020



Use distilled water  
Używaj wody destylowanej  
Verwenden Sie destilliertes Wasser  
Utilizzare acqua distillata  
Utilizar agua destilada  
Folositi apă distilată  
Utiliser de l'eau distillée  
Gebruik gedestilleerd water  
Используйте дистиллированную воду

 **Bestel hier!**

 **PRAXISDIENST**  
Medische artikelen sinds 1953

Manufacturer: Enbio Group AG, Eichengasse 3, 4702 Oensingen, Switzerland







**NEW**  
**NOWY**  
**NEU**  
**NUOVO**  
**NUEVO**  
**NOU**  
**NOUVEAU**  
**NIEUW**  
**НОВЫЙ**



**HALF USED**  
**ZUŻYTY W POŁOWIE**  
**HALB GEBRAUCHT**  
**USATO A METÀ**  
**MEDIO USADO**  
**JUMATATE FOLOSITĂ**  
**MOITIÉ UTILISÉE**  
**HALF VERBRUIKT**  
**ПОЛОВИНА**  
**ИСПОЛЬЗОВАНИЯ**



**BUY NEW**  
**KUP NOWY**  
**NEU KAUFEN**  
**COMPRARE NUOVO**  
**COMPRAR NUEVO**  
**CUMPĂRĂ NOU**  
**ACHETER NOUVEAU**  
**KOOP NIEUW**  
**КУПИТЬ НОВЫЙ**



 **PRAXISDIENST**  
Medische artikelen sinds 1953

 **Bestel hier!**