



## IH 58 Nebuliser


**german  
engineering**


Nebulisation of liquid medicine using compressed air technology



Small, light and silent

For the treatment of the upper and lower respiratory tract

For use in colds, asthma, respiratory disease, etc.

- for fast and effective inhalation

High proportion of respirable particles

High nebulisation capacity (approx. 0.25 ml/min.)

- short inhalation time

Accessories:

Nebuliser bottle, mouthpiece, nosepiece, adult mask, child mask, compressed air tube, mains adapter with micro-USB cable, filter, storage bag and carrying bag

Can be disinfected

Operating pressure: approx. 0.25-0.5 bar

Medical device

Mains operation:

Input: 100-240 V~; 50-60 Hz; 0.5 A

Output: 5 V  $\equiv$  2 A

Product measurements (LxWxH): 110 x 62 x 47 mm

Product weight: approx. 176 g

5 year guarantee

Sales unit: 3 / Shipping carton: 4 x 3

EAN no.: 4211125 60202 4

Item no.: 602.02



**PRAXISDIENST**  
 Medical Supplies since 1953  
