

TEST TUBE WITHOUT RIM



DESCRIPTION

This test tube is made of high-quality borosilicate glass. The high material durability to heat, sudden temperature changes and most chemicals makes this test tube ideal for a variety of laboratory procedures. These test tubes are available in different sizes.

PRODUCT DETAILS

- Test tube without rim
- Material: 3.3 borosilicate
- Colour: transparent
- Temperature resistant
- Resistant to most chemicals
- Can be autoclaved
- Diameter (depending on size): 12 mm, 15 mm or 16 mm
- Height (depending on size): 75 mm, 100 mm or 150 mm
- Quantity: 1 box with 100 pieces

Overview

Article Number:
T135842-T135844



Category:
Lab Equipment

 **PRAXISDIENST**
Medical Supplies since 1953

 **Order here!**

Teqler is a trademark of NetMed S.à.r.l.

 NetMed S.à.r.l.
8, Am Scheerleck, 6868 Wecker, Luxembourg

 +352 267149 09  +352 267149 19  info@netmed.lu

Social Media

