

## TEQLER CHILDREN MASK WITH MOTIVE



### DESCRIPTION

The 3-layer Teqler children's mouthguard with motif has elastic ear loops and an adjustable nose clip. This ensures that the mask is very comfortable to wear and can be optimally adapted to the child's face. The children's mask complies with DIN EN 14683 type IIR and has a bacteria filtering capacity of 99%. The mouthguard is optionally available with animal or airplane motif.

### PRODUCT DETAILS

- Teqler children's mouthguard with motif
- 3-ply
- With elastic ear loops and adjustable nose clip
- Bacteria filtering capacity of 99
- Conform to DIN EN 14683 TYPE IIR
- Comfortable and tight fit
- CE certified medical device
- Available with either airplane or animal motifs
- For single use
- Pack of 50 masks

### Overview



Article Number:  
**T143476 - T143477**

Category:  
**Surgical Masks**



Teqler is a trademark of NetMed S.à.r.l.

 NetMed S.à.r.l.  
8, Am Scheerleck, 6868 Wecker, Luxembourg

 +352 267149 09  +352 267149 19  [info@netmed.lu](mailto:info@netmed.lu)

Social Media

