

DERMLITE DL5 TECHNICAL DATA AND COMPARISON TABLE

	DL1	DL100	DL200 HR	DL200 Hybrid	DL3N	DL4W	DL4	DL5
POLARIZATION	Polarized & Non-Pol.	Polarized	Polarized	Polarized & Non-Pol.	Polarized & Non-Pol.	Polarized & Non-Pol.	Polarized & Non-Pol.	Variable polarization (incl. polarized/non-polarized toggle)
OPTICS	ø15mm 10x	ø15mm 10x	ø25mm 10x	ø25mm 10x	ø25mm 10x	ø30mm 10x	ø30mm 10x	ø32mm 10x
ILLUMINATION	4 white LEDs (pol./non-pol.)	8 white pol. LEDs	21 white pol. LEDs	21 white LEDs (15 pol., 6 non-pol.)	18 white LEDs (12 pol., 6 non-pol.), 10 orange pol. LEDs	24 white LEDs (18 pol., 6 non-pol.)	18 white LEDs (12 pol., 6 non-pol.), 6 orange pol. LEDs	16x white, 4x UV (365nm), 4x Pigment Boost, 1 torch LED
INFECTION CONTROL	-	-	IceCap®	IceCap	-	IceCap	IceCap	IceCap
POWER	Rechargeable	Disposable (Type 2CR5)	Rechargeable	Rechargeable	Rechargeable, includes charging base	Rechargeable	Rechargeable	Rechargeable, includes charging base
SMARTPHONE USE	Yes	-	Yes	Yes	Yes	Yes	Yes	Yes, includes MCC adapter
CONTACT PLATE	Bayonet	-	Snap-on	Snap-on	Fixed	Snap-on	Snap-on	Bayonet style

INCLUDED ACCESSORIES	MCC smartphone adapter		Silicone Sleeve with Lanyard, IceCap supply, USB cable	Silicone Sleeve with Lanyard, IceCap supply, USB cable	Desktop Charger Leather Pouch Eyepiece	USB Cable Belt Holster Silicone Sleeve IceCap sample supply	USB Cable Belt Holster Silicone Sleeve IceCap sample supply	Desktop charging base, phone adapter, magnetic eyepiece, 100 mm ruler, leather belt pouch, protective holster, 5 IceCaps, 2m USB-A to USB-C charging cable, microfiber cleaning cloth
----------------------	--	--	--	--	--	--	--	---



PRAXISDIENST
Medical Supplies since 1953



Order here!