



Veterinary HCT Meter

Calibrated and validated specifically
for use in dogs, cats and horses



**QUICK, EASY,
HCT RESULTS.
2 YEAR WARRANTY.**

"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

WOODLEY
EQUIPMENT COMPANY LTD.

InSight™ Veterinary HCT Meter

The InSight HCT meter is a low cost handheld meter for the measurement of Canine, Feline and Equine haematocrit and total haemoglobin in capillary and venous whole blood.

KEY FEATURES

- Calibrated for Dogs, Cats and Horses
- Accurate and precise
- Handheld
- 30 second test time
- One drop whole blood (Venous, Capillary or EDTA)
- Stores 250 results
- HCT results reported as %
- Room temperature strip storage



Calibrated and validated specifically for use in dogs, cats and horses



KIT CONTENTS

- InSight HCT Meter
- InSight HCT Test Strips (10)
- InSight HCT Code Cards - Cat, Dog & Horse
- Lancets (10)
- Lancing Device
- Instruction Manual
- CR2032 Batteries (2)
- Carry Case

SPECIFICATIONS

- **Parameters:** Haematocrit and Haemoglobin
- **Measuring Principle:** Optical reflectance
- **Strip:** InSight HCT Test Strips
- **Display:** LCD
- **Weight:** <58g (Batteries included)
- **Dimensions:** 102 x 50 x 19mm
- **Power:** DC 6V (Two CR2032 Batteries)
- **Operating Temperature:** 15-30°C
- **Reportable Range:** 12-72%

+44 (0) 1204 669033 (Option 1)

sales@woodleyequipment.com

www.woodleyequipment.com

Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

**"PASSIONATE ABOUT
LABORATORY DIAGNOSTICS"**

**WOODLEY
EQUIPMENT COMPANY LTD.**