



ENCORE[®] Latex Micro

S U R G I C A L G L O V E S

Latex, powder-free surgical glove with exceptional dexterity and tactile sensitivity for surgical finesse

KEY FEATURES AND BENEFITS:

- 20% thinner on average than Ansell's standard latex powder-free surgical gloves with soft latex formulation for unsurpassed dexterity and tactile sensitivity
- Its brown color helps reduce glare of surgical lights and makes it suitable to be used as an underglove
- The DERMASHIELD™ Technology internal coating facilitates dry or damp donning and intra-operative changes

CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.

Recommended For

- All surgeries, especially delicate procedures in ophthalmology, micro, and cardiovascular
- Damp donning
- Enhanced tactile sensitivity and dexterity, especially for delicate procedures
- Prevention of eye fatigue

FEATURED TECHNOLOGIES



ENCORE[®] The Confidence to Perform

Dedicated to delivering optimized solutions engineered to provide a high degree of confidence to surgeons, O.R. nurses and surgical technicians seeking protection with a high level of dexterity, comfort and grip so they can perform at their best.

ENCORE[®] Latex Micro

SURGICAL GLOVES

PRODUCT INFORMATION

Material	Natural Rubber Latex
Powder Content	Powder-Free
Color	Brown
Thickness	Micro
Cuff Style	Beaded with SUREFIT™ Technology
External Glove Surface	Smooth with Micro-Textured Finish
Internal Glove Surface	Polymer Coated with DERMASHIELD™ Technology
Allergy Prevention	None
Grip Level	Moderate
Double Gloving Recommendation	Outer glove or underglove
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)
Freedom from Holes (Inspection level I)	0.65 AQL
Product Standards	ASTM D3577, ASTM D7160, EN 16523-1, EN 374:2003, EN 420:2003 + A1:2009, EN 455 1-4, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	Zinc Dibutyldithiocarbamate (ZDBC), Zinc Diethyldithiocarbamate (ZDEC)

ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	330104055	330104060	330104065	330104070	330104075	330104080	330104085	330104090

Contact your Ansell representative for ordering or more information.

North America
 US Tel: 800 952 9916
 CA Tel: 844 494 7854
 Email: insidesalesus@ansell.com
 Sandel US Tel: 866 764 3327
 Sandel email: sandel-cs@ansell.com

Central & South America
 Tel: +52(442) 296 20 50
 Email: cslac@ansell.com
Brazil
 Tel: +55-11-3356-3100
 Email: luvass.medicas@ansell.com

Europe, Middle East & Africa
 Tel: +32 (0) 2 528 74 00
 Email: info@ansell.eu

Asia Pacific
 Tel: +603 8310 6688
 Email: apac.medical@ansell.com

India
 Tel: +91 22 38124500
 Email: info.india@ansell.com
Australia & New Zealand
 Tel: +61 3 9270 7270
 Email: protection@ap.ansell.com
China
 Tel: +86(21) 3827 5005
 Email: infochina@ansell.com